



DELIVERING
INNOVATION
TO PROTECT

 extract-technology.com

AFTERSALES SOLUTIONS



Extract Technology have been working with the worlds leading pharmaceutical and chemical companies for over 35 years so you can be sure we have the knowledge, expertise and resources to help you continue to provide a safe working environment.

That is why we place a huge emphasis on the maintenance and upkeep of your equipment. We can achieve this through the supply of approved spare parts, regular servicing and the facility to upgrade and relocate your equipment whenever necessary.



Extract Technology's Multi Port Glove Test Kit


The failure of component parts can cause costly disruption to production schedules, so it is essential to have quick access to replacement parts.

SPARE & REPLACEMENT PARTS

AT EXTRACT TECHNOLOGY WE STOCK AN EXTENSIVE RANGE OF SPARE PARTS AND CONSUMABLES AVAILABLE TO DISPATCH ANYWHERE IN THE WORLD AT SHORT NOTICE.

OUR RANGE OF REPLACEMENT PARTS INCLUDES:

-> Replacement filters
-> Cooling coils and control valves
-> Safe change filter bags and rings
-> Gloves (Hypalon, Butyl, Neoprene)
-> Ceiling and bleed PLF screens
-> Gaskets
-> Glove/Sleeve combinations
-> Replacement glazing panels
-> Beta transfers (Sleeves, containers)
-> Magnehelic gauges





OUR RANGE OF VALIDATION SERVICES INCLUDES:

- ...> Downflow Booths
- ...> Aseptic / Flexible Isolators
- ...> Fume Cabinets/Ductless Fume Cabinets
- ...> Laminar Flow Benches/ Enclosures
- ...> LEV Systems
- ...> Clean Rooms
- ...> Biological Safety Cabinets Class I / Class II/ Class III
- ...> Positive & Negative Pressure Systems
- ...> Forensic Downflow Tables
- ...> P2 & 3 Laboratories

PLANNED PREVENTATIVE MAINTENANCE

OUR PLANNED PREVENTATIVE MAINTENANCE PACKAGE IS DESIGNED TO HELP SUSTAIN THE OPTIMUM PERFORMANCE OF YOUR EQUIPMENT BOTH NOW AND INTO THE FUTURE.

IN TODAY'S ENVIRONMENT OF STRICT REGULATION AND COMPLIANCE, THE NEED TO GUARANTEE THAT YOUR EQUIPMENT IS PERFORMING CORRECTLY IS CRITICAL. EXTRACT TECHNOLOGY CAN HELP YOU ACHIEVE THIS WITH REGULAR, SCHEDULED EXAMINATION AND TESTING



PLANNED PREVENTATIVE MAINTENANCE PACKAGE



Extract's Planned Preventative Maintenance package provides a scheduled visit to undertake tests and inspection of equipment, along with necessary minor repairs and adjustments. This ensures equipment is maintained and disruption to production schedules through failure is minimised.

All tests carried out to current ISO14644 standard and aligned to the HSG 258 (Controlling airborne contaminants at work 2017) recorded on quality controlled documentation with approved test methods.

- ...> All tests performed in accordance with industry standard testing requirements and recorded on quality control documents
- ...> Certificate of Conformity issued for Health and Safety records
- ...> Planned, coordinated visits to ensure minimum disruption
- ...> All examination, maintenance and testing carried out by trained professionals using calibrated instruments
- ...> Clear reports on findings
- ...> Honest assessment of equipment with recommendation on upgrades

EQUIPMENT UPGRADES AND RELOCATIONS

EXTRACT TECHNOLOGY'S HIGHLY SKILLED INSTALLATION TEAMS CAN HELP WITH ANY RELOCATION REQUIREMENTS, ENSURING EQUIPMENT IS MOVED SAFELY AND WITH MINIMUM DISRUPTION.

Additionally, upgrading your containment equipment is possible with a number of options including: replacement wall panels, containment screens, safe change filtration and cooling packages.

OUR RANGE OF UPGRADE AND RELOCATIONS SERVICES INCLUDES:

- > Equipment Relocations
- > Cooling Coil Packages
- > Solutions to Reduce Operator Exposure Levels
- > Facility & Airlock Modifications
- > Versatile 5 Dimensional Rigid Containment Screens
- > Safe Change & Gel Seal Filter Upgrades
- > Multi Ports Glove Test Kits

Other Services

OPERATOR TRAINING/ RE-TRAINING

Extract Technology offer a full and effective training and re-training programme for both operators and maintenance personnel in the safe use of the equipment and SOP's.

PERSONAL AIR SAMPLING/ PARTICLE COUNTING

We can undertake Operator Exposure Monitoring using either real time aerosol monitoring or personal sampling. Additionally, airborne particle count testing can be carried out to determine area classification.

FILTER INTEGRITY TESTING

Independent HEPA filter testing can be carried out to the current legislation ensuring that all installations are documented as compliant.

RE-VALIDATION IQ/OQ

A full re-validation package is available, from initial document creation to full on-site execution. These range from a basic service to full installation and operational qualification (IQ & OQ) to suit the customer's needs.

DOCUMENTATION

Our Compliance department can create any type of documentation from simple pre-delivery documents to full system qualifications. Operation & Maintenance manuals include all specification sheets and engineering drawings.



WHERE WE WORK

- Europe
- Asia
- Australasia
- Russia
- Middle East
- Africa
- South America
- Central America
- North America

FOR MORE INFORMATION ON OUR PRODUCTS PLEASE VISIT OUR WEBSITE OR CONTACT US DIRECTLY:

EXTRACT TECHNOLOGY

 extract-technology.com
 info@extract-technology.com



EXTRACT TECHNOLOGY LTD UK, EUROPEAN AND ROW OFFICE

SPARES - +44 (0) 1484-317-033
SERVICE - +44 (0) 1484-354-499
UPGRADES - +44 (0) 1484-354-493

BRADLEY JUNCTION IND ESTATE
LEEDS RD, HUDDERSFIELD
HD2 1UR
UNITED KINGDOM

EXTRACT TECHNOLOGY AMERICAS OFFICE

PARTS - + 1 608-562-7749
SERVICE - + 1 608-562-7711

618 STATE STREET
NEW LISBON, WI 53950
U.S.A.