



extract
TECHNOLOGY



REPLACEMENT PARTS | SERVICING | UPGRADES

Aftermarket Support




Extract Technology have been working with the worlds leading pharmaceutical and chemical companies for over 30 years so you can be sure we have the knowledge, expertise and resources to help you continue to provide a safe working environment.

That is why we place a huge emphasis on the maintenance and upkeep of your equipment. We can achieve this through the supply of approved spare parts, regular servicing and the facility to upgrade and relocate your equipment whenever necessary.




FOR MORE INFORMATION ON AFTERMARKET SUPPORT CALL
+44 (0) 1484 432 727 OR EMAIL service@extract-technology.com



extract
TECHNOLOGY

Spares & Replacement Parts



All types of filters supplied



Monitoring gauges & consumables



Perfect laminar flow screens



Gloves, gauntlets & gloveports

The failure of component parts can cause costly disruption to production schedules, so it is essential to have quick access to replacement parts.

At Extract Technology we stock an extensive range of spare parts and consumables available to dispatch anywhere in the world at short notice.

Our range of replacement parts include:

- Replacement filters
- Cooling coils and control valves
- Safe change filter bags and rings
- Gloves (Hypalon, Butyl, Neoprene)
- Ceiling and bleed PLF screens
- Gaskets
- Glove/Sleeve combinations
- Replacement glazing panels
- Beta transfers (Sleeves, containers)
- Magnehelic gauges



Safe change bags and rings



Cooling Coil

Servicing & Remedial Work



Planned Preventative Maintenance

Our Planned Preventative Maintenance package is designed to help sustain the optimum performance of your equipment both now and into the future.

In today's environment of strict regulation and compliance, the need to guarantee that your equipment is performing correctly is critical.

Extract Technology can help you achieve this with regular, scheduled examination and testing.



Our range of validation servicing includes:

- Downflow Booths
- Aseptic / Flexible Isolators
- Fume Cabinets/ Ductless Fume Cabinets
- Laminar Flow Benches/ Enclosures
- LEV Systems
- Clean Rooms
- Biological Safety Cabinets Class I / Class II/ Class III
- Positive & Negative Pressure Systems
- Forensic Downflow Tables
- P2 & 3 Laboratories



Downflow Booth

Extract's Planned Preventative Maintenance package provides a scheduled visit to undertake tests and inspection of equipment, along with necessary minor repairs and adjustments. This ensures equipment is maintained and disruption to production schedules through failure is minimised.

- All tests performed in accordance with industry standard testing requirements and recorded on quality control documents
- Certificate of Conformity issued for Health and Safety records
- Planned, coordinated visits to ensure minimum disruption
- All examination, maintenance and testing carried out by trained professionals using calibrated instruments
- Clear reports on findings
- Honest assessment of equipment with recommendation on upgrades

Equipment Upgrades & Relocations

Extract Technology's highly skilled installation teams can help with any relocation requirements, ensuring equipment is moved safely and with minimum disruption.

Additionally, upgrading your containment equipment is possible with a number of options including, replacement wall panels, containment screens, safe change filtration and cooling packages etc.

Our upgrade and relocations packages includes:

- Equipment Relocations
- Cooling Coil Packages
- Solutions to Reduce Operator Exposure Levels
- Versatile 5 Dimensional Rigid Containment Screens
- Safe Change & Gel Seal Filter Upgrades
- Facility & Airlock Modifications



Equipment relocations



Cooling control



5D Containment Screens

Other Services

OPERATOR TRAINING/ RE-TRAINING

Extract Technology offer a full and effective training and re-training programme for both operators and maintenance personnel in the safe use of the equipment and SOP's.

PERSONAL AIR SAMPLING/ PARTICLE COUNTING

We can undertake Operator Exposure Monitoring using either real time aerosol monitoring or personal sampling. Additionally, airborne particle count testing can be carried out to determine area classification.

FILTER INTEGRITY TESTING

Independent HEPA filter testing can be carried out to the current legislation ensuring that all installations are documented as compliant.

RE-VALIDATION IQ/OQ

A full re-validation package is available, from initial document creation to full on-site execution. These range from a basic service to full installation and operational qualification (IQ & OQ) to suit the customer's needs.

DOCUMENTATION

Our Compliance department can create any type of documentation from simple pre-delivery documents to full system qualifications. Operation and Maintenance manuals include all specification sheets and engineering drawings.

EXTRACT TECHNOLOGY LIMITED

Bradley Junction Industrial Estate, Leeds Road,
Huddersfield, HD2 1UR, United Kingdom

T. +44 (0) 1484 432 727
F. +44 (0) 1484 432 659
W. www.extract-technology.com

